




NEAR-INFRA RED SPECTROMETER

Part Number: NIR256-1.7T1

- 900 - 1700nm
- InGaAs Array
- T.E. Cooled Array
- Portable For Lab or Field Use
- No Moving Parts
- 100% Permanent Wavelength Calibration*



CDI NIR spectrometer use T.E. cooled, high performance Linear InGaAs arrays to provide complete spectra in milliseconds.

CDI spectrometer are ideal for process, QC, laboratory and educational purposes.

We also offer intergrated systems in portable computers, industrial chassis's with single board compters and battery powered portables.

SEMICONDUCTORS - TELECOMMUNICATIONS - POLYMER - THIN FILMS
PHARMACEUTICAL - AGRICULTURAL - PETROCHEMICAL - CHEMICAL

*Winthin .5 nm with AutoCal option



Variable Path Transmission Dip Probe



2633 Foundation Drive, South Bend, IN 46628
(574) 288-7338 Fax (574) 288-7339
Email sales@controldevelopment.com
www.controldevelopment.com

Control Development, Inc. also manufactures integrated spectroscopy systems for: UV-VIS, NIR, RAMAN and OEM customers.



NIR



PORTABLE



RAMAN



MY INSTRUMENT



OEM